



## NanoTack® Flex Anatomic Labrum Restoration Systems

A New Trajectory in Hip Labrum Restoration

### Product Highlights

- 25° delivery system may facilitate increased safety margins by reducing joint penetration via a flexible delivery system allowing multiple angles of approach
- Compact 1.4mm, high-strength implant allows implantation near the labro-cartilage junction
- 2-stage fixation yields comparable strength to larger hip anchors\*
- Anchor and pilot hole trajectory are matched allowing for seamless implantation where intended

\*Pivot Medical Internal Data PR01569-A Section 4.1.1 and PR01618-A Section 4.6.1.1



*NanoTack Flex Anatomic Labrum Restoration System.*

## NanoTack Flex Anatomic Labrum Restoration System

### Ordering Information

Description	Catalog No.
NanoTack Suture Anchor, 1.4mm with Flex Inserter	CAT01858
NanoTack Flex Drill Bit	CAT01859
NanoTack Flex Guide	CAT01860
NanoTack Flex Obturator	CAT01861
NanoTack Guide	CAT01371
NanoTack Obturator	CAT01372
KnotPusher	CAT01518

**For more information about the NanoTack Flex System,  
or to request a demo, please call Pivot Medical at  
1-866-596-2022**



*The NanoTack Flex System allows placement near the labro-cartilage junction in multiple trajectories through a straight or curved guide allowing placement of multiple non-everting stitches to allow for complete labrum refixation even in hard to reach areas.*

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: NanoTack, Stryker, The Movement in Hip Restoration. All other trademarks are trademarks of their respective owners or holders.

**Pivot Medical, Inc.**  
A Stryker Sports Medicine Company  
247 Humboldt Court  
Sunnyvale, CA 94089 USA